



Medical Imaging Resources, Inc.

120 Enterprise Drive

Ann Arbor, Michigan 48103

Phone: 888.323.1316

Fax: 734.426.2003

www.mobileleasing.com

www.medimatingsales.com



MAGNETOM Espree – the first Open Bore MR We see a way to provide patients with CT-like comfort in a 1.5T MRI 70 cm + 125 cm + 1.5T and Tim – a combination never seen before in MRI. Which is why we hesitate to even call it that. In truth, it's an entirely different breed. MAGNETOM Espree combines Open Bore MRI with Tim (Total imaging matrix) technology. **MAGNETOM Espree – the second Tim system**

It combines Open Bore MRI with Tim (Total imaging matrix) technology. [Click here to see why Tim has changed MR forever. see more>>](#) **MAGNETOM Espree – the first Open Bore MR**

The challenge wasn't to make a more "open" Open. The challenge was to retain the performance of a 1.5T system in a combination with the most open design possible. This combination of 1.5T, Open Bore, and Tim is something only Siemens could do. An innovation unlike anything on the market today. [See more >>](#) **MAGNETOM Espree – Revolutionary workflow**

The most intelligent system through integrated workflow automation. [See more >>](#) **System variants of MAGNETOM Espree**

MAGNETOM Espree, Tim [32 x 8] or MAGNETOM Espree, Tim [76 x 18] For studies ranging from clinical routine up to advanced applications. It utilizes highest SNR for superb image quality and comes with a Z-engine (33 mT/m, SR 100T/m/s) gradient system.

All attempts have been made to ensure accurate data. Medical Imaging Resources, Inc. assumes no responsibility for any unintentional errors or omissions.