



**Medical Imaging Resources, Inc.**

120 Enterprise Drive

Ann Arbor, Michigan 48103

Phone: 888.323.1316

Fax: 734.426.2003

---

[www.mobileleasing.com](http://www.mobileleasing.com)

[www.medimagingsales.com](http://www.medimagingsales.com)



## Upgrade MAGNETOM Symphony with **Tim**

Join the Revolution. More than just a routine hardware a quantum leap to leading edge MRI technology and MAGNETOM Symphony and Tim – a perfect combination patient friendliness and fast scanning capabilities.

Tim, Total imaging matrix, is a completely new, clinically proven technology. This Tim upgrade provides a unique opportunity to integrate the latest technology into your MAGNETOM Symphony:

- A complete new RF chain with up to 32 RF channels, the best technology in the industry for high speed Parallel Imaging.
- A set of Matrix coils can flexibly combine up to 76 coil elements for unlimited Parallel Imaging.
- The new telescopic patient table for new applications, from extended Field of View imaging up to whole-body examinations.

All attempts have been made to ensure accurate data. Medical Imaging Resources, Inc. assumes no responsibility for any unintentional errors or omissions.