



Medical Imaging Resources, Inc.
120 Enterprise Drive
Ann Arbor, Michigan 48103
Phone: 888.323.1316
Fax: 734.426.2003

www.mobileleasing.com
www.medimagingsales.com



MAGNETOM Trio MAGNETOM Trio

MAGNETOM Trio MAGNETOM Trio is the first generation of true 3T whole-body 3T MR system. The largest diversity of iPAT-compatible coils combined with the standard 8 RF channels ensure the best 3T examinations for all organs.

The guaranteed homogeneity on 40 cm FoV, including in the z-direction and the standard high order shim enable best spine and abdominal imaging ever in the coronal plane.

MAGNETOM Trio opens new doors to applications in abdominal, cardiac, spine and orthopaedic MR at 3T

All attempts have been made to ensure accurate data. Medical Imaging Resources, Inc. assumes no responsibility for any unintentional errors or omissions.